

Senator Mould Chart

senator *the difference is in the smile*



UPPER ANTERIORS

Square Moulds



Ovoid Moulds



Tapering Moulds



LOWER ANTERIORS

3C

34.0
5.0
8.0



3N

33.0
4.5
9.0



2D

33.0
5.0
9.5



2P

33.5
5.0
9.5



2E

33.0
5.0
9.0



47

36.5
5.5
9.5



3P

38.0
5.5
10.5




26

35.5
5.0
11.0




4H

36.5
5.0
9.0



3R

40.0
5.5
12.0



Senator conforms to BS EN ISO: 22112:2006



Mould 4H
Width of 6 on flat 46.5
Width of left central incisor 8.5
Length of left central incisor 10.5



Mould 4H
Width of 6 on flat 36.5
Width of left central incisor 5.0
Length of left central incisor 9.0



Mould 29M
Width of 4 on flat U - 29.5
Depth of crown - left first molar 7.5
U signifies maxillary dentition



Mould 29M
Width of 4 on flat L - 31.5
Depth of crown - left first molar 8.5
L signifies mandibular dentition

All measurements in millimetres

Senator Articulation Table

Anterior		Posteriors				
U	L	Functional		20	33	SB
3C	3C	1	11S/13S	29M/L	16	
2D	2D	1	11S/13S	29M/L	16	
3D	3C	1	11S/13S	29M/L	16	
5P	3N	2	11S/13S	29M/L	18	18SB
2E	2E	2	11S/13S	31M/L	18	18SB
3X	4H	4	11S/13S	31M/L	20	20SB
1N	3N	4	11S/13S	31M/L	18	18SB
2P	2P	4	11S/13S	31M/L	18	18SB
A25	2E	2	11S/13S	31M/L	18	18SB
A26	47	4	11S/13S	31M/L	18	18SB
4H	4H	4	11S/13S	31M/L	18	18SB
136	4H	4	11S/13S	31M/L	20	20SB
137	47	4	11S/13S	31M/L	20	20SB
3N	3N	4	11S/13S	31M/L	18	18SB
3P	3P	4	11S/13S	31M/L	20	20SB
266	26	4	11S/13S	31M/L	20	20SB
267	3R	4	11S/13S	31M/L	20	20SB

POSTERIORIORS

Functional Moulds

33° Moulds

1
U 26.5
8.0
L 28.5
7.75



16
U 29.0
8.25
L 32.0
9.0



2
U 29.0
9.0
L 32.0
8.0



18
U 31.5
9.0
L 35.0
10.0



4
U 30.0
9.0
L 33.0
7.0



20
U 33.5
9.5
L 36.5
9.5



11S
U 29.5
8.0
L 31.0
7.0



18SB
U 31.5
9.0
L 35.0
10.0



13S
U 31.0
8.5
L 32.5
7.25



20SB
U 33.5
9.5
L 36.5
9.5



20° Moulds

29M

U 29.5
7.5
L 31.5
8.5



29L

U 29.5
8.25
L 31.5
8.5



31M

U 31.0
8.0
L 34.0
8.5



31L

U 31.5
8.5
L 34.0
9.0



Zero Degree Moulds

29Z

U 29.0
7.0
L 31.0
8.0

